beurer

BF 1000 Super Precision

DE	Diagnosewaage	
	Gebrauchsanweisung	2
EN	Diagnostic bathroom scale	
	Instructions for use.....	20
FR	Pèse-personne impédancemètre	
	Mode d'emploi.....	38
ES	Báscula diagnóstica	
	Instrucciones de uso	57
IT	Bilancia diagnostica	
	Istruzioni per l'uso	75
TR	Diyanoz terazisi	
	Kullanım kılavuzu.....	93
RU	Диагностические весы	
	Инструкция по применению	110
PL	Waga diagnostyczna	
	Instrukcja obsługi	129

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Read these instructions for use carefully. Observe the warnings and safety notes. Keep these instructions for use for future reference. Make the instructions for use accessible to other users. If the device is passed on, provide the instructions for use to the next user as well.

Contents

1. Included in delivery	21	9. Analysis of measured values	28
2. Signs and symbols	21	10. Evaluating the results	31
3. Warnings and safety notes	22	11. Transferring measurements via app	33
4. Device description	24	12. Additional information	34
5. Initial use	25	13. Cleaning and maintenance	34
6. Configuring the diagnostic bathroom scale with an app	26	14. Disposal	35
7. Alternatively: configuring the diagnostic bathroom scale without an app	26	15. What if there are problems?	35
8. Carrying out the measurement	27	16. Technical specifications	36
		17. Warranty / service	37

Getting to know your device

Functions of the device

The purpose of this digital scale is weighing and calculating your personal fitness data. It is intended for private use.

The scale is equipped with the following functions, which can be used by up to 10 different people:

- Body weight measurement
- Body fat percentage calculation
- Visceral fat percentage calculation (organ fat/abdominal fat)
- Body water content calculation
- Muscle percentage calculation
- Bone mass calculation
- Basal metabolic rate and active metabolic rate calculation.

The scale also includes the following additional features:

- Cross measurement with 8 electrodes allows for segment analysis, i.e. a separate display of individual body parts (arms, legs, torso)
- Innovative menu navigation on large display (in 5 languages),
- Switch between kilograms "kg", pounds "lb" and stones "st"
- Battery replacement display for low batteries
- Automatic user recognition
- Storage of the last 30 measurements for 10 users if the data cannot be transferred to the app
- Guest mode (without memory function).

The measurement principle

This scale works on the principle of B.I.A. (bioelectrical impedance analysis). This involves the calculation of body content in a matter of seconds, using a current that cannot be felt, is completely harmless and does not pose any risk.

When this measurement of the electrical resistance (impedance) is considered alongside constants and/or individual values (age, height, gender, activity level), the body fat percentage and other variables in the body can be calculated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as the fat cells and bones are barely able to conduct the current due to very high levels of resistance.

The electrical contact to the body is enhanced through 4 hand electrodes and 4 foot electrodes to achieve a higher measurement accuracy and better analysis. This sensitive analysis technology means readings may change slightly when measurements are repeated. This is caused by changes during the impedance measurement (e.g. different electrical contact to the electrodes by the hands and feet, different distribution of water in the body).

Please note that the values calculated by the diagnostic bathroom scale only represent an approximation of the actual medical analytical values for the body. Only a specialist doctor can precisely establish the body fat, body water, muscle percentage and bone structure using medical methods (e.g. computer tomography).

1. INCLUDED IN DELIVERY

Check that the exterior of the cardboard delivery packaging is intact and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Service address.

- Diagnostic bathroom scale with raised display on aluminium stand
- 4 x 1.5V type AA batteries
- These instructions for use
- 2 x quick guide
- Mounting kit (4 screws with 2 replacement screws + Allen key)

2. SIGNS AND SYMBOLS

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

	Warning Warning instruction indicating a risk of injury or damage to health
	Important Safety note indicating possible damage to the device/accessory
	Product information Note on important information
	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.
	CE labelling This product satisfies the requirements of the applicable European and national directives.
	The products demonstrably meet the requirements of the Technical Regulations of the EAEU
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Do not dispose of batteries containing hazardous substances with household waste
	United Kingdom Conformity Assessed Mark

	Importer symbol		
	Manufacturer		Separate the product and packaging elements and dispose of them in accordance with local regulations.
	Separate the packaging elements and dispose of them in accordance with local regulations.		Marking to identify the packaging material. A = Material code, B = Material number: 1-7 = Plastics, 20-22 = Paper and cardboard
	Consult instructions for use		Danger of slipping! Do not step on the scale if the surface or your feet are wet
	Place the scale on a stable, level surface. Do not place the scale on a carpet.		Danger of tipping! Do not stand on the outer edge.

3. WARNINGS AND SAFETY NOTES

WARNING

- The device is only intended for weighing humans and for recording your personal fitness data. The device is only intended for private use and not for medical or commercial purposes. The device must only be used for the purpose for which it was designed and in the manner specified in these instructions for use.
- **The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.**
- **Do not use the device during pregnancy.**
- **Danger of tipping!** Never place the scale on an uneven surface. Never step onto the outer edge of the scale on just one side.
- **Danger of slipping!** Do not step onto the scale with wet feet. Do not step onto the scale when the surface is wet. Do not step onto the scale wearing socks.
- Please note that measuring tolerances are possible for technical reasons as this scale is not calibrated for use in a professional medical context.
- Obtain professional advice from appropriately qualified personnel (e.g. a doctor) when undertaking any treatment or diet to address being overweight or underweight. The health values determined using the scale can assist the qualified personnel in their diagnoses.
- Inaccurate or implausible results may occur when determining body fat in children, performance athletes or people with medical/physical limitations (e.g. oedema symptoms, osteoporosis, people who are taking cardiovascular medication or medication which changes the vascular system, people with major anatomical abnormalities on the legs).
- Repairs may only be carried out by Beurer Customer Services or authorised retailers. Never open the scale yourself.

Safety notes on handling batteries

- Always insert the batteries correctly, taking into account the polarity (+ / -). Keep batteries clean and dry and away from water. Always select the correct battery type.
- Never short-circuit batteries and battery compartment contacts.
- Never charge, forcibly discharge, heat, disassemble, deform, encapsulate or modify batteries.
- Never weld or solder on batteries.

- Never mix batteries of different manufacturers, capacities (new and used), size and type within a device.
- **Risk of explosion!** Failure to comply with the points mentioned above can result in personal injury, overheating, leakage, venting, breakage, explosion or fire.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- **Choking hazard!** Keep batteries out of the reach of children. Seek medical attention immediately if swallowed.
- Never allow children to replace batteries without adult supervision.
- Store batteries away from metal objects in a well-ventilated, dry and cool place.
- Never expose batteries to direct sunlight or rain.
- Remove batteries from the device if it is not going to be used for a long period of time.
- Dispose of discharged batteries immediately and properly. Never dispose of batteries in fire.
- When disposing of batteries, keep batteries with different electrochemical systems separate.


General notes

- The scale's maximum capacity is 200 kg (440 lb/31 st). The results for weight measurement are displayed in increments of 50 g (0.1 lb).
- Measurements of body fat, body water and muscle content are displayed in increments of 0.1%.
- The calorie requirement is indicated in increments of 1 kcal.
- When supplied to the customer, the scale is set to weigh and measure in "cm" and "kg". You can change the unit settings during initial use with the app.

General tips


- As far as possible, weigh yourself at the same time every day (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be carried out using measurement with the hands and bare feet; the soles of the feet may be lightly dampened where appropriate.
- Results may be unsatisfactory if the soles of the feet are completely dry or have a large amount of hard skin since conductivity will be impaired.
- Stand up straight and still during measurement. Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in your body can disperse.
- Remember that only the long-term trend is important. Short-term changes in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.

4. DEVICE DESCRIPTION

1. Display
2. Menu: right button (sensor button)
3. Electrodes, right hand
4. Menu: down button (sensor button)
5. Menu: up button (sensor button)
6. Set button, confirm the entry (sensor button)
7. ON/OFF button (sensor button)
8. Electrodes, left hand
9. Menu: left button (sensor button)
10. Aluminium stand
11. Foot electrodes, right foot
12. Foot electrodes, left foot

Display


1. Date	2. Time
3. <i>Bluetooth®</i>	4. General overview
5. Segment analysis	6. Graphic representation
7. User memory	8. Average values
9. Calorie consumption	10. Body fat percentage with interpretation
11. Body water content with interpretation	12. Muscle percentage with interpretation
13. Bone mass with interpretation	14. Body weight with BMI interpretation
15. Weight unit (kg, lb or st)	

5. INITIAL USE


Positioning the scale


Place the scale on a flat and stable surface. A solid surface is a prerequisite for a precise measurement.


1. Open the cardboard box.


2. Remove the polystyrene cover on the left and right.


3. Lift up the upper and lower part of the aluminium stand at the same time.


4. Insert the upper part of the aluminium stand into the lower part and secure with screws at the top.


5. Position the aluminium stand upright in the cardboard box. Slide the scale to the edge...

6. ...and secure with screws at the bottom.

7. Hold the stand in place with one hand and remove the film with the other.

8. Insert batteries the correct way round in the underside of the scale (+/-). Stand the scale upright on a flat surface.

6. CONFIGURING THE DIAGNOSTIC BATHROOM SCALE WITH AN APP

To calculate your body fat percentage and additional body values, you must save the personal user data on the scale.

The scale has 10 user memory spaces which allows, for example, you and your family members to save personal settings.

It is also possible to activate users via other mobile end devices on which the “beurer HealthManager Pro” app has been installed or by changing the user profile in the app (see the app settings).

Keep the smartphone close to the scale in order to maintain an active *Bluetooth®* connection during initial use.

- **Activate *Bluetooth®* in your smartphone's settings.**
- **Download the free "beurer HealthManager Pro" app from the Apple App Store (iOS) or from Google Play (Android™).**
- **Start the app and follow the instructions.**
 - The following settings must be set or entered in the “beurer HealthManager Pro” app:

User data	Set values												
Initials	Max. 3-digit (A-Z, 0-9)												
Gender	Male (♂), female (♀)												
Height	100 to 250 cm (3' 03" to 8' 02")												
Date of birth	Year, month, day												
Activity level	The medium-term and long-term perspective is decisive when selecting the activity level:												
	<table><tr><td></td><td>Physical activity</td></tr><tr><td>1</td><td>None.</td></tr><tr><td>2</td><td>Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).</td></tr><tr><td>3</td><td>Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.</td></tr><tr><td>4</td><td>High: physical exertion, at least 4 to 6 times a week for 30 minutes each.</td></tr><tr><td>5</td><td>Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.</td></tr></table>		Physical activity	1	None.	2	Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).	3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.	4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.	5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.
		Physical activity											
	1	None.											
	2	Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).											
	3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.											
4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.												
5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.												
Target weight	kg/lb/st (for progress graph)												

- The app assigns the next free user memory space (person 1 to person 10) to you.
- Carrying out the measurement: For automatic personal recognition, the first measurement must be assigned to your personal user data. The app will prompt you to do this. For how to do this, please see chapter “8. Carrying out the measurement”.

7. ALTERNATIVELY: CONFIGURING THE DIAGNOSTIC BATHROOM SCALE WITHOUT AN APP

We recommend using the app for the initial use and implementing all settings. See chapter “6. Configuring the diagnostic bathroom scale with an app”. However, this can also be carried out on the device itself.

Specify device/time/user settings on the control unit. To do so, use the [⏮] and [⏭] buttons to make your selection and confirm with [SET] respectively.

a.) Setting the time

[SET] → “SET TIME” → [SET]

Sequence	Set values
Date	Year–Month–Day
Time	Hours–Minutes

b.) Specifying device settings

[SET] → “SET UNIT” → [SET]

Sequence	Set values
Language	DE EN FR ES IT
Measurement unit	kg–lb–st
Time	24 h–12 h
Guest mode	ON–OFF
Bluetooth®	ON–OFF  When set to “OFF”, a <i>Bluetooth</i> ® connection is only possible if the scale is switched on. This can be useful for extending the battery life or to prevent undesirable radio radiation.

c.) Setting user data

[SET] → “SET USER” → [SET]

User data	Set values
Initials	Max. 3-digit (A-Z, 0-9)
Gender	Male () , female ()
Height	100 to 250 cm (3' -03" to 8' -02")
Date of birth	Year, month, day
Activity level	1 to 5 (see section 6)
Target weight	kg/lb/st (for progress graph)
User PIN	Assigned automatically for <i>Bluetooth</i> ® use at a later time. This information is only required for app set-up in the future, if existing scale users are to be added during app installations at a later time.

Logging in to the app at a later time

If you have carried out the initial use of the scale without using the app, you can add your existing scale user to the “beurer HealthManager Pro” app. This requires you to enter your 4-digit user “PIN”. This process is described in the app.

8. CARRYING OUT THE MEASUREMENT

8.1 Body weight measurement

- Press the  button on the control unit. „0.0“ kg appears in the display.
- The scale is then ready to take a weight measurement. To measure your weight, step onto the scale with your weight distributed evenly on both feet. The weight data will not be stored; user allocation will not be possible, as there is no BIA measurement via the hands or feet.


8.2 Carrying out guest measurement

For one-time measurement or if no user data is to be stored. For this, guest mode must be switched on in the device settings.

[] → “Guest” → [SET] or switch scale on

User data	Set values
Gender	Male () , female ()
Height	100 to 250 cm (3' -03" to 8' -02")
Date of birth	Year, month, day
Activity level	1 to 5 (see section 6)

8.3 User measurement, total body analysis

- Press the [] button on the control unit. „0.0“ kg appears in the display.  If guest mode has been activated in the scale settings, step onto the scale when it is switched off. - Stand on the scale barefoot holding the control unit in your hands. Hold the control unit with the underside of your thumbs on the front electrodes and with your fingers pressed firmly against the back electrodes. - Stand still with your weight distributed equally on both feet. Make sure there is no skin contact between your two feet, legs, calves and thighs.	
- If the scale recognises you automatically, the impedance for the analysis will be measured after the weight measurement. This is indicated by the symbol showing current passing through the body. - As soon as the total body analysis is complete, the overview will be shown. The measurement is automatically stored.	

Note:

The measurement will not be saved if the following applies:

- The user is not automatically detected. This happens if your measured weight differs from the previous measurement by +/- 3 kg or if another user has similar values (+/- 3 kg) to yours. In this case, press the [] button to select the user manually and to carry out the measurement again.
- The user is standing on the scale with shoes on.
- The user is not holding the handheld unit with their hands.

9. ANALYSIS OF MEASURED VALUES


The results will be analysed automatically after taking the measurement or after selecting the user manually on the control unit. Select your user using the [] button. Confirm using the [SET] button or wait 3 seconds.

General overview


The general overview appears at a glance in the display either before or after a total body analysis for your user.

Segment analysis for body fat and muscle percentage


a) Body fat in left and right arm

- Press the menu button [] until is selected. The segment analysis for body fat in the left and right arm appears in the display.


b) Body fat in left and right leg

- Use the [] button to select the detailed display of body fat for the left and right leg.


c) Body fat in the torso

- Use the [] button to select the detailed display of body fat in the torso.


d) Visceral fat (organ fat/abdominal fat)

- Use the [] button to select the detailed display of visceral fat.


a) Muscle percentage in left and right arm

- Use the [] button to select the detailed display of muscle percentage in the left and right arm.


b) Muscle percentage in left and right leg

- Use the [] button to select the detailed display of muscle percentage in the left and right leg.

g) Muscle percentage in torso

- Use the [] button to select the detailed display of muscle percentage in the torso.


No visceral fat or muscle percentages in specific segments of the body (arm/leg/torso) will be displayed for users aged 10 – 17 years.

Graphic representation

a) Body weight

- Press the menu button [>] until is selected. The trend for your body weight appears in the display. The trend for your body weight is shown from measurement to measurement. Later measurements are added on the right. The reference line shows your target weight.

You can change your target weight in the user data at any time.

b) Body fat

- Use the [<] button to switch to showing body fat values .

The trend for your body fat values is shown from measurement to measurement. The reference line represents the average value for all body fat measurements; the bars show whether the values have gone up or down.

c) Water proportion

- Use the [<] button to switch to showing water proportion .

The trend for your body's water proportion is shown from measurement to measurement. The reference line represents the average value for all water proportion measurements; the bars show whether the values have gone up or down.

d) Muscle percentage


- Use the [<] button to switch to showing muscle percentage .

The trend for your muscle percentage is shown from measurement to measurement. The reference line represents the average value for all muscle percentage measurements; the bars show whether the values have gone up or down.

M Showing stored measurements

- Press the menu button [>] until is selected. The most recent measured stored measurement appears in the display. The measurement date and the memory space number are shown at the top, e.g. (03-03) means the third memory space = the last of a total of 3 occupied memory spaces.
- Use the [<], [>] buttons to cycle through the memory spaces.


Average values

- Press the menu button [>] until is selected. All average values in the stored measurements are shown in the display.

The dates of the first and the last measurement are shown at the top.

Body water content

The body water content in % is normally within the following ranges:

Male				Female			
							
Age	Poor	Good	Very good	Age	Poor	Good	Very good
10–100	<50 %	50–65 %	>65 %	10–100	<45 %	45–60 %	>60 %

Body fat contains relatively little water, therefore body water content may be below the standard values in persons with a high body fat percentage. On the other hand, endurance athletes may exceed the standard values due to having low levels of fat and a high muscle percentage. The body water calculation performed using this scale is not suitable for drawing medical conclusions, for example in the case of age-related water retention. Consult your doctor where necessary. A high body water content is generally desirable.

Muscle percentage

The muscle percentage measurement measures the skeletal muscle mass. No organ muscles are taken into account. The muscle percentage is normally within the following ranges:

Male				Female			
							
Age	Low	Normal	High	Age	Low	Normal	High
10-14	<44%	44–57 %	>57 %	10-14	<36%	36–43 %	>43 %
15-19	<43%	43–56 %	>56 %	15-19	<35%	35–41 %	>41 %
20-29	<42%	42–54 %	>54 %	20-29	<34%	34–39 %	>39 %
30-39	<41%	41–52 %	>52 %	30-39	<33%	33–38 %	>38 %
40-49	<40%	40–50 %	>50 %	40-49	<31%	31–36 %	>36 %
50-59	<39%	39–48 %	>48 %	50-59	<29%	29–34 %	>34 %
60-69	<38%	38–47 %	>47 %	60-69	<28%	28–33 %	>33 %
70-100	<37%	37–46 %	>46 %	70-100	<27%	27–32 %	>32 %

Bone mass

Like the rest of our body, our bones are subject to natural growth, shrinking and ageing processes. Bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass then begins to reduce. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeleton can be further increased by means of targeted muscle build-up. Note that this scale identifies the mineral content of the bone. It is very difficult to influence the bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

Important

Please do not confuse bone mass with bone density. The bone density can only be determined by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions (e.g. lying in bed 24 hours a day). This value is primarily dependent on weight, height and age. It is displayed on the diagnostic scale in the unit kcal/day and calculated using the scientifically recognised Mifflin St. Jeor equation.

This is the amount of energy that is required by your body under all circumstances and must be re-supplied to the body in the form of food. If your energy intake is below this level in the longer term, this can result in damage to health.

AMR

The active metabolic rate (AMR) is the amount of energy that the active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the diagnostic bathroom scale calculates this by means of the activity level (1–5) that has been entered.

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink. If less energy is taken in than is used over a long period of time, the body primarily takes the difference from the fat stores with resulting weight loss. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

Temporal context of the results

Please note that only the long-term trend is significant. Brief deviations in weight within a few days are normally caused by loss of fluids.

The interpretation of the results is based on changes in total body weight, percentage of body fat, body water and muscle content, as well as on the length of time over which these changes occur. Rapid changes within the scope of a few days are to be considered separately from medium-term changes (in the scope of weeks) and long-term changes (months).

It can be said as a basic rule that short-term changes in weight are almost entirely changes in water content, whereas medium-term and long-term changes may also relate to the fat percentage and the muscle percentage.

- If your weight drops in the short term but your body fat percentage increases or stays the same, you have only lost water, for example following a training session, visit to the sauna or a diet aimed only at fast weight loss.
- On the other hand, if your weight increases in the medium term and your body fat percentage drops or stays the same, you may have built up valuable muscle mass.
- If your weight and body fat percentage both fall at the same time, then your diet is working—you are losing fat mass.
- Ideally, you should support your diet with physical activity, fitness or strength training. This enables you to increase your muscle percentage in the medium term.
- Body fat, body water and muscle percentages should not be totalled (muscle tissue also contains components made of body water).

11. TRANSFERRING MEASUREMENTS VIA APP

If the app is not open, the newly assigned measurements are saved on the scale. A maximum of 30 measurements per user can be saved on the scale. The saved measurements are transferred automatically to the app when you open the app within the *Bluetooth®* range. It is possible to automatically transfer data when the scale is switched off if *Bluetooth®* is enabled in the device settings.

12. ADDITIONAL INFORMATION

Deleting scale data

a) Delete user with all settings and saved values

- After switching on, select your user by pressing the [] button. If necessary, press the button several times and confirm with [SET].
- When **0.0** appears in the display, hold down the [SET] button for 12 seconds to delete the user.
- Use the [, ] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete the user, confirm “NO” using the [SET] button.

All saved values and settings for this user will be deleted.

b) Delete individual saved value for a user

- Select a measured value in the user memory [] and briefly press the [SET] button on the control unit to delete the stored value displayed.
- You can also do this immediately after the measurement is taken.
- Use the [, ] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete the value, confirm “NO” using the [SET] button.

c) Delete all stored values for a user

- Open the user memory [] and press the [SET] button on the control unit for around 5 seconds to delete all values.
- Use the [, ] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete, confirm “NO” using the [SET] button.

Changing the batteries

The control unit has a low battery indicator. “LOW BATTERY” appears in the display if the batteries in the scale are too weak. In this case, the batteries should be replaced. Connect the scale to the app again. The system time is updated.

Note:

- When replacing the batteries, use batteries of the same type, make and capacity.
- Do not use rechargeable batteries.
- Use batteries free from heavy metals.

13. CLEANING AND MAINTENANCE

Clean the surface of the scale from time to time using a slightly damp cloth. Do not use the scale again until the surface is completely dry.

IMPORTANT

- Do not use any abrasive cleaning products.
- Never submerge the scale in water. Never rinse the scale under running water.
- Do not place any objects on the scale when it is not in use.
- Protect the scale from knocks, moisture, dust, chemicals, drastic changes in temperature and nearby sources of heat (ovens, heaters).

14. DISPOSAL

Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

For environmental reasons, do not dispose of the device in the household waste at the end of its service life. Dispose of the device at a suitable local collection or recycling point in your country.

Dispose of the device in accordance with EC Directive–WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

15. WHAT IF THERE ARE PROBLEMS?

If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Solution
„ERROR:_OVERLOAD“	Maximum weight capacity of 200 kg exceeded.	Loads must be limited to 200 kg.
„LOW_BATTERY“	Batteries are nearly flat.	Replace the batteries (see section on changing the batteries).
„MEMORY_FULL“	30 measurements have been stored.	Transfer the measured values.
„ERROR:_CONTACT“ „ERROR“	The electrical resistance between the electrodes and the soles of the feet is too high (e.g. with very hard skin).	Repeat the measurement barefoot. Dampen the soles of your feet slightly. If necessary, remove hard skin on the soles. Hold the hand electrodes with the fingers behind and the balls of the thumb in front.
„ERROR:_DATA“	The BIA measurement was successful, but the fat percentage is outside the measurable range (less than 3% or greater than 65%).	Please repeat the measurement barefoot, or with the soles of the feet lightly dampened if necessary. If necessary, change your set activity level.
Empty display	Batteries are completely flat.	Replace the batteries (see section on changing the batteries).
	The batteries are not inserted properly.	Check for correct polarity, (see “Initial use” -> “Changing the batteries”).
Only weight is displayed, no diagnosis	Unknown measurement as it is outside the user assignment limit or a unique assignment is not possible.	Select user on the scale using the [] button or repeat user assignment in the app. Hold the electrodes of the handheld unit and stand barefoot on the scales.

Incorrect weight measurement

Possible faults	Solution
Incorrect zero point of scale.	Switch the scale on, wait approx. 5 seconds and repeat the measurement.
Scale is on a carpet.	Place the scale on a flat, solid surface. Stand on the outside of the scale as far as possible without it tipping.

Possible faults	Solution
Incorrect weight is displayed.	Stand as still as possible.
Incorrect weight is displayed.	Wait until the scale switches itself off. Then repeat the measurement.

No **Bluetooth®** connection (📶 symbol is missing)

Possible faults	Solution
Radio connection	The range in the open is approx. 10 m. Walls and ceilings reduce the range. Other radio waves may interfere with the transmission. Therefore, do not position the scale near devices such as WIFI routers, micro-waves, induction hobs etc.
Bluetooth® connection	Completely close the app (including in the background). Switch Bluetooth® off and back on again. Switch the smartphone off and back on again. Remove the batteries in the scale for a short period and reinsert them. See FAQs at www.beurer.com .
Connection set-up taking too long	Faster Bluetooth® interval when the scale is switched on.

16. TECHNICAL SPECIFICATIONS

Scale dimensions	380 x 390 x 1000 mm
Scale batteries	4 x 1.5V type AA
Measurement range	8 kg to 200 kg
Scale interval d:	The display can be read in increments of 0.05 kg.
Repetition accuracy:	The measuring tolerance for repeated measurements is +/- 0.3 kg (several measurements in direct succession on the same scale with position of scale and person the same as far as possible).
Absolute precision:	In comparison with a calibrated weight, the measured value is +/-0.5% +0.05 kg. E.g. at 50 kg this corresponds to +/- 0.3 kg; at 100 kg this corresponds to +/-0.55 kg.
Data transfer via Bluetooth® wireless technology	This product uses Bluetooth® low energy technology, frequency band 2.402–2.480 GHz, transmission power max. 0 dBm, compatible with Bluetooth® 4.0 smartphones/tablets

System requirements and list of compatible devices for the “beurer HealthManager Pro” app*

Subject to technical changes.

*This product satisfies the requirements of the applicable European directives.
We hereby confirm that this product complies with the European RED Directive 2014/53/EU.
The CE Declaration of Conformity for this product can be found under:
<https://www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php>

17. WARRANTY / SERVICE

Further information on the warranty and warranty conditions can be found in the warranty leaflet supplied.

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The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.

Beurer GmbH • Söflinger Str. 218 • 89077 Ulm (Germany) • www.beurer.com
www.beurer-gesundheitsratgeber.com • www.beurer-healthguide.com

Beurer UK Ltd., Suite 16, Stonecross Place, Yew Tree Way
WA3 2SH Golborne, United Kingdom